




OSS-DBS: Quick Start

By Konstantin Butenko kbutenko@bwh.harvard.edu

Why OSS-DBS

- Advanced E-field modeling:
 - tissue dispersiveness and capacitance
 - conductivity tensors
 - adaptive mesh refinement
- Cable equations
- Supports pathway activation modeling
- For the theory, refer to <https://doi.org/10.1371/journal.pcbi.1008023>
- Example of use <https://doi.org/10.1016/j.nicl.2022.103185>



Why not:

- Computationally demanding (do not run on your old laptop)
- Advanced simulations require understanding of model parameters


— Overview

- Installation
- Setting up in Lead-DBS
 - *VTA approximation*
 - *Pathway activation modeling (PAM)*
 - *Multiple protocols*
- OSS-DBS GUI
- OSS-DBS terminal


Installation

OSS-DBS is distributed as an external library in Lead-DBS

And runs within *Docker* (can be thought as a virtual environment)



- Install *Docker*
 - Windows and macOS users should install *Docker Desktop* and provide corresponding permissions
 - Linux users are referred to <https://docs.docker.com/desktop/install/linux-install/>, make sure you install the appropriate version
- Enable development environment in *Preferences File*
- In MATLAB command line, run `ea_checkOSSDBSInstall`
 - Follow prompted suggestions, you might need to update *pip*



Installation

Without Lead-DBS installation routines (for confident Linux users)

- Install *Docker* following <https://docs.docker.com/desktop/install/linux-install/>
- Install/update *pip*, run `python3 -m pip install PyQt5` and `python3 -m pip install h5py`
 - I had to re-install some PyQt5 dependencies on *Ubuntu 22.04*
- Create a docker group, add your user and reboot
 - `sudo groupadd docker`
 - `sudo usermod -aG docker <USERNAME>`
- Pull the docker image `sudo docker pull ningfei/oss-dbs:latest`
- Open a terminal in `leaddbs/.../ext_libs/OSS-DBS/` and run the following to allow non sudo execution of docker container
 - `docker build --build-arg UID=$(id -u) --build-arg GID=$(id -g) -t ningfei/oss-dbs:custom -f custom.Dockerfile .`
- Run OSS-DBS with the following command
 - `docker run --name OSS_docker --volume <OSS-DBS path>:/opt/OSS-DBS --volume <patient stim folder path>:/opt/Patient --cap-add=SYS_PTRACE -it --rm <docker image>`
 - E.g. for a test patient *pBK*: `docker run --name OSS_docker --volume ~/Documents/leaddbs/ext_libs/OSS-DBS:/opt/OSS-DBS --volume ~/Documents/Data/pBK:/opt/Patient --cap-add=SYS_PTRACE -it --rm ningfei/oss-dbs:custom`

Additional files

- Lead-DBS *data* folder (*templates/...*) contains tensor data for MNI space (*IITMeanTensor.nii.gz* & *IITMeanTensor_NormMapping.nii.gz*)

- Advance users can use custom DTI data

- Fiber atlases are stored as '*lead-dbs-folder/connectomes/dMRI/'atlas_name'/data.mat*


- Pathway atlases are stored as

'lead-dbs-folder/connectomes/dMRI_MultiTract/'connectome_name/'pathway_name.mat'

- Example of a fiber atlas:

https://github.com/SFB-ELAINE/OSS-DBS/blob/67a30ac26b300e7517b5c302912d7130c126db45/OSS_platform/Example_fibers.mat

- Contact us for connectomes, we have plenty!



Setting up in Lead-DBS

Required minimum for non-experts

- Template or native space (the latter is preferred)
- Source settings either in Lead-DBS or OSS-DBS (the latter is easier)
- Use *settings* to specify the simulation setup

Stimulation Parameters

Patient: CLOVER05

Stimulation: 20201220182053 Settings

Model: OSS-DBS (Butenko 2020)

Beta ☒ Estimate in template space ☐ Add StimSet

Right Hemisphere

Sources	3	mA	0	V
0	V	0	V	

	%		%
K3	0	K7	50
K2	0	K6	0
K1	50	K5	0
K0	0	K4	0
Case	100		

OSS-DBS setting

Settings for OSS-DBS (Butenko et al., 2020)

☐ Calculate Axon Activation

Multi-Tract: Petersen_HDP_GPI_Th

Axon Length: 10 mm

Fiber Diameter: 5.7 mm

E-Field Threshold Presets: Petersen_HDP_GPI_Th

E-Field Threshold: 0.323 V/mm

☒ Use tensor data: IITmean_tensor.nii.gz

☒ Interactive mode

Save



Check if OSS-DBS
GUI will be used

External DTI data

Setting up in Lead-DBS: VTA approximation

Standard approach of VTA approximation by $|E|$ -threshold is supported

- VTA is seeded with MRI resolution
- Maximum $|E|$ during the pulse
- Rattay's function upon request
- Directional VTA in development




Setting up in Lead-DBS: PAM

Activation of ‘realistic’ axons

Choose a connectome/pathway

If one pathway, settings here

- The same morphology within one pathway
- When modeling one pathway, store a file with fibers as
‘lead-dbs-folder’/connectomes/dMRI/‘pathway_name’/data.mat
- When modeling multiple pathways, store them as
.../dMRI_MultiTract/‘connectome_name’/‘pathway_name.mat’
 - **IMPORTANT:** Review ‘Setting’ when modeling multiple pathways!
- Fibers are pre-filtered using intentionally exaggerated Kuncel VTA



Setting up in Lead-DBS: PAM

Activation of ‘realistic’ axons

- Axons are seeded around active contacts
- Axon length: limit for long fibers to reduce computational costs
- Fiber diameter: educated guess
- Axon model: McIntyre’s mammalian axon (or classic McNeal’s model upon request)

Choose a connectome/pathway

If one pathway, settings here



OSS-DBS setting

Settings for OSS-DBS (Butenko et al., 2020)

☒ Calculate Axon Activation

Multi-Tract: Petersen_HDP_GPI_Th

Axon Length: 10 mm

Fiber Diameter: 5.7 mm

E-Field Threshold Presets:

E-Field Threshold: 0.323 V/mm

☒ Use tensor data: IITmean_tensor.nii.gz

☒ Interactive mode

Save

Setting for Connectome

	Axon Length (mm)	Fiber Diameter (mm)
M1_hdp_face	10.0000	5.7000
M1_hdp_lowerex	10.0000	5.7000
M1_hdp_upperex	10.0000	5.7000
PreMotor_hdp	10.0000	5.7000
SMA_hdp_lhs	10.0000	5.7000
ansa_lenticularis	15.5000	11.5000
lenticular_fasciculus	15.5000	11.5000


Save

Setting up in Lead-DBS: PAM

Activation of 'realistic' axons

- Purple - damaged, red - activated (i.e. AP in response to the DBS pulse)
- Damaged due to the intersection with the electrode, encapsulation and CSF
- Activation is also stored in Network_state.h5 file (dataset for each projection)
- Fiber activation is displayed automatically in electrode-scene

Show fiber/axon status →




← Show fibers /axons

Today	Today	Today
20201220182053	<ul style="list-style-type: none"> Allocated_axo...N_nodes.csv Allocated_axons.h5 Animation_Field_in_time Axons_in_time Boston_Scien...la_position.py Brain_substitute.brep CSF_ref fiberActivatio...hdp_face.mat fiberActivatio...p_lowerex.mat fiberActivatio...p_upperex.mat fiberActivatio...otor_hdp.mat Field_solutions Field_solutions_functions Images Meshes MRI_DTL_derived_data Neuron_model_arrays oss-dbs_parameters.mat Results_adaptive Results_rh 	<ul style="list-style-type: none"> Activation_All...ularis_338.csv Activation_All...culus_380.csv Activation_All...p_face_19.csv Activation_All...werex_55.csv Activation_All...perex_22.csv Activation_All...r_hdp_27.csv Axon_state_a...nticularis.mat Axon_state_l...fasciculus.mat Axon_state_...hdp_face.mat Axon_state_...p_lowerex.mat Axon_state_...p_upperex.mat Axon_state_P...otor_hdp.mat Connection_s...nticularis.npy Connection_s...fasciculus.npy Connection_s...hdp_face.npy Connection_s...p_lowerex.npy Connection_s...upperex.npy Connection_s...otor_hdp.npy Last_run_in_...lenticularis.csv Last_run_in_l...fasciculus.csv

Setting up in Lead-DBS: PAM

Activation of 'realistic' axons

- Purple - damaged, red - activated (i.e. AP in response to the DBS pulse)
- Damaged due to the intersection with the electrode, encapsulation and CSF
- Activation is also stored in Network_state.h5 file (dataset for each projection)
- Fiber activation is displayed automatically



Setting up in Lead-DBS: Multiple protocols

Field superposition instead of recalculations

$$\nabla \cdot \left[\left(\sigma(\mathbf{r}, \omega) + j\omega\epsilon(\mathbf{r}, \omega) \right) \nabla \phi \right] = 0$$

- Linearity of Laplace's equation allows to 'combine' and scale solutions (el. potential distribution) for each active contact
- Straight-forward for constant current mode
- StimSet input in mA(!)
- Use case grounding


	Contact0	Contact1	Contact2	Contact3	Contact4	Contact5	Contact6	Contact7
1	-1.2000	-1.2000	1.5000	NaN	1.5000	-1.2000	-1.2000	NaN
2	2.1000	NaN	-2.1000	NaN	2.1000	NaN	-2.1000	NaN
3	1.3000	NaN	1.3000	1.3000	NaN	1.3000	NaN	NaN
4	NaN	NaN	2.2000	NaN	-2.2000	NaN	2.2000	NaN

	Contact0	Contact1	Contact2	Contact3	Contact4	Contact5	Contact6	Contact7
1	3.5000	NaN	NaN	-1.2000	NaN	NaN	-1.2000	1.5000
2	NaN	1.7000	NaN	NaN	NaN	-1.7000	NaN	NaN
3	1.5000	1.5000	1.5000	NaN	NaN	-4.5000	NaN	NaN

Setting up in Lead-DBS: Multiple protocols


Field superposition instead of recalculations

- Linearity of Laplace's equation allows to 'combine' and scale solutions (el. potential distribution) for each active contact
- Straight-forward for constant current mode
- StimSet input in mA(!)
- Use case grounding




OSS-DBS GUI

Allows to access more modeling parameters



Lead-DBS input →



OSS-DBS GUI

Deeper understanding of the model

OSS-DBS v0.3 for Lead-DBS

Simulation State
 Simulation Setup
 Visualization

 Advanced Simulation Parameters

MRI and DTI Data

Geometries

Coordinates from ea_reconstruction.mat

Electrode Type	Boston_Scientific_Vercise_Cartesia	?
Implantation Coordinate X	11.0679	mm ?
Implantation Coordinate Y	-13.2106	mm ?
Implantation Coordinate Z	-9.6288	mm ?
2nd point on lead X	12.5273	mm ?
2nd point on lead Y	-10.2542	mm ?
2nd point on lead Z	-3.7323	mm ?
Turn around lead's axis	7.57	deg ?

Current-Controlled Mode ☒ ?

Signal Vector [-0.003, None, None] ?

Solver Type GMRES ?

OSS-DBS v0.3

Elaine
electrically active implants
SIB 1270

Icons by [Adioma](#)

Load
 Save As
 Reset
 Run

OSS-DBS v0.3 for Lead-DBS

Simulation State
 Simulation Setup
 Visualization

 Advanced Simulation Parameters

MRI Data ready ☐ ?

DTI Data ready ☐ ?

Initial Neuron Array ready ☐ ?

Geometry and Initial Mesh ready ☐ ?

Adjusted Neuron Array ready ☐ ?

Signal Generation done ☐ ?

CSF Refinement done ☐ ?

Adaptive Refinement done ☐ ?

Computations in Spectrum done ☐ ?

Continue Interrupted Computations ☐ ?

Scaling and IFFT done ☐ ?

Skipping steps if interrupted

OSS-DBS v0.3


Elaine
electrically active implants
SIB 1270

Icons by [Adioma](#)




Load
 Save As
 Reset
 Run

OSS-DBS GUI


Deeper understanding of the model



Parameters of the electrode-tissue interface



Detailed description of the DBS pulse




OSS-DBS GUI

Deeper understanding of the model

- Input configuration is stored as a python dictionary `GUI_inp_dict.py` in *stim* folder of the patient


Load a preconfigured simulation



Resave your configuration as this file to use in the simulation or save in another .py file for future

Instead of using GUI, you can modify entries of the dictionary directly (for advanced users)

Launch a simulation in a Docker container



OSS-DBS outputs

Behind the scene

- Various outputs and intermediate files are stored in patient folders
- Can be rather large (depending on the neuron array)
- What is interesting: *Tensors/*, *Images/*, *Field_solutions/*
- Status files: *success_rh*, *skip_rh*, *fail_rh* (same with *_lh*)

OSS-DBS terminal

Behind the scenes

- Shows the simulation flow (but concisely)
- Closes automatically if everything is fine
- Log files are available in stim. Folders use them to report errors

Check out the simulation progress in complete_log...

```
Setting up the simulation...
Processing |#####| 3/11
Preparing neuron models and FEM mesh...
/opt/OSS-DBS/OSS_platform
Processing |#####| 4/11
Generating geometry and mesh...
Processing |#####| 5/11
Adjusting neuron models
Processing |#####| 6/11
Generating DBS pulse...
```

```
[2022-04-14 19:20:00,858]:Processed frequencies:
[2022-04-14 19:20:00,858]:0.0 130.0 221.92388155425118
[2022-04-14 19:20:34,656]:313.84776310850236 497.69552621700467 865.3910524340093
[2022-04-14 19:21:09,605]:1600.7821048680187 3071.5642097360374 6013.128419472075
[2022-04-14 19:21:45,272]:11896.25683894415 23662.5136778883 47195.0273557766
[2022-04-14 19:22:18,300]:94260.0547115532 188390.1094231064 376650.2188462128
[2022-04-14 19:22:18,627]:--- Sorting the obtained solution
[2022-04-14 19:22:19,626]:Saved sorted solution in Field_solutions/
[2022-04-14 19:22:19,626]:----- Parallel calculations took 4 min 4 sec -----

[2022-04-14 19:22:19,626]:----- Conducting signal scaling and IFFT -----
[2022-04-14 19:22:55,980]:25% of neuron models were processed
[2022-04-14 19:23:27,081]:51% of neuron models were processed
[2022-04-14 19:24:01,817]:75% of neuron models were processed
[2022-04-14 19:24:33,769]:----- Signal scaling and IFFT took 2 min 14 sec -----

[2022-04-14 19:24:33,769]:----- Calculating impedance -----

[2022-04-14 19:24:33,792]:Max impedance: 1531.155924348973
[2022-04-14 19:24:36,867]:----- Estimating neuron activity -----
[2022-04-14 19:25:15,613]:0 models were activated
[2022-04-14 19:25:15,616]:0.0% activation (including damaged neurons)

[2022-04-14 19:25:15,622]:----- NEURON calculations took 0 min 38 sec -----

[2022-04-14 19:25:15,622]:----- Simulation run took 9 min 13 sec -----
```

OSS-DBS support

- Do not hesitate to contact us (kbutenko@bwh.harvard.edu)
- The software is rather new, not all features are implemented
- Cluster solution is available upon request